



# WHY YOU CAN'T GET A PLANAR FLAP FROM THE HANSATOME® OR MORIA® M2



DIFFERING SPEEDS IN THE ARC  
CREATE UNEVEN FLAP THICKNESS



CONSTANT LINEAR SPEED  
CREATES A PLANAR FLAP

The Hansatome® microkeratome is a product and registered trademark of Bausch & Lomb, Inc. The M2 microkeratome is a product of and Moria® is a registered trademark of Moria Surgical, S.A. MED-LOGICS, Inc. is not affiliated with Bausch & Lomb, Inc. or Moria Surgical, S.A.



**Your patients deserve better.**

