

unique  **edge**[®]
MICROSURGICAL KNIVES



Precision you can see. Quality you can feel.



Clinicians know that minimizing risks and providing positive patient outcomes are the most important issues to consider when choosing the right surgical instruments. Unique Technologies offers a line of UniqueEdge[®] knives that consistently provide both the quality and precision necessary to achieve these goals.

PRECISION YOU CAN SEE

- Our vacuum-purified double melt, hardened and tempered stainless steel technology produces a uniform strength and ductility that is critical to creating our precise cutting edge for smooth, uniform incisions
- Consistent blade dimensions
- Flawless cutting edge
- Superior wound architecture
- Minimal wound leakage

QUALITY YOU CAN FEEL

- Effortless tissue penetration
- Excellent tactile response with minimal drag
- Consistent reproducible incisions

PERFORMANCE YOU CAN TRUST



- Accurate, well controlled incisions
- Tightly sealing incisions
- Sutureless closure
- Rapid wound healing
- Reduced potential for infection

TEXTURED HANDLE

- Market proven design
- Enhanced instrument control
- Excellent tactile response

MATTE FINISHED BLADE

- Improves view of the operating field
- Eliminates reflective glare from operating room light sources
- Maximal visibility of the knife within the eye chamber






PRECISION YOU CAN SEE

The microscope reveals the flawless, finished cutting edge that makes our microsurgical blades instruments of precision.

BROAD CUSTOMIZATION CAPABILITIES

With the wide variety of blade dimensions and geometries available, the UniqueEdge® line is ready to meet your most demanding surgical requirements.

You can operate with complete confidence with a knife that flawlessly compliments your procedure, equipment, instruments, and surgical technique.



Stab Knives



ITEM NO.	WIDTH	DESCRIPTION
52-1501	15°	Full Handle
52-2201	22.5°	Full Handle
52-3001	30°	Full Handle
52-4501	45°	Full Handle
52-1531	15° 3mm	Full Handle, restricted depth blade
52-1551	15° 5mm	Full Handle, restricted depth blade
52-1532	15° 3mm	Blade and tip, restricted depth blade
52-1552	15° 5mm	Blade and tip, restricted depth blade
52-1502	15°	Blade and tip
52-2202	22.5°	Blade and tip
52-2241	22.5° 4mm	Full Handle, restricted depth blade
52-3036	30° 3.5mm	Full Handle, restricted depth blade

Clear Corneal Knives



ITEM NO.	WIDTH	DESCRIPTION
52-2563	2.5mm	Angled, full handle
52-2663	2.65mm	Angled, full handle
52-2763	2.75mm	Angled, full handle
52-2863	2.85mm	Angled, full handle
52-3063	3.0mm	Angled, full handle
52-3263	3.2mm	Angled, full handle

Clear Corneal Knives- Side Cutting/ Parallel Sides



ITEM NO.	WIDTH	DESCRIPTION
52-2635	2.65mm	Angled, full handle
52-2735	2.75mm	Angled, full handle
52-2835	2.85mm	Angled, full handle
52-3035	3.0mm	Angled, full handle
52-3235	3.2mm	Angled, full handle

Phaco Slit Knives



ITEM NO.	WIDTH	DESCRIPTION
52-2561	2.5mm	Angled, full handle, bevel up
52-2661	2.65mm	Angled, full handle, bevel up
52-2761	2.75mm	Angled, full handle, bevel up
52-2861	2.85mm	Angled, full handle, bevel up
52-2961	3.0mm	Angled, full handle, bevel up
52-3261	3.2mm	Angled, full handle, bevel up
52-3561	3.5mm	Angled, full handle, bevel up
52-5261	5.2mm	Angled, full handle, bevel up
52-2901	3.0mm	Straight, full handle, single bevel
52-3201	3.2mm	Straight, full handle, single bevel
52-2762	2.75mm	Angled, blade and tip, bevel up
52-2862	2.85mm	Angled, blade and tip, bevel up
52-2962	3.0mm	Angled, blade and tip, bevel up
52-3262	3.2mm	Angled, blade and tip, bevel up
52-2831	2.85mm	Angled, full handle, double bevel
52-2931	3.0mm	Angled, full handle, double bevel

MicroPhaco Slit Knives



ITEM NO.	WIDTH	DESCRIPTION
52-1463	1.4mm	Angled, full handle, bevel up
52-1563	1.5mm	Angled, full handle, bevel up
52-1863	1.8mm	Angled, full handle, bevel up
52-2063	2.0mm	Angled, full handle, bevel up
52-2261	2.2mm	Angled, full handle, bevel up
52-2263	2.2mm	Angled, full handle, double bevel

Please call your UniqueEdge® distributor to inquire about knife designs or sizes not listed.

Precision you can see. Quality you can feel.

Crescent Knives



ITEM NO.	WIDTH	DESCRIPTION
54-1060	1.25mm	Angled, full handle, bevel up
54-1010	2.0mm	Angled, full handle, bevel up
54-1000	2.0mm	Straight, full handle
54-5062	2.0mm	Angled, blade and tip, bevel up
54-5000	2.0mm	Straight, blade and tip
54-1013	2.3mm	Angled, full handle, bevel up

Vitrectomy Knives



ITEM NO.	WIDTH	DESCRIPTION
51-1901	19 gauge	Straight, full handle
51-2001	20 gauge	Straight, full handle
51-2061	20 gauge	Angled, full handle, bevel up
51-2301	23 gauge	Straight, full handle
51-2361	23 gauge	Angled, full handle, bevel up
51-2501	25 gauge	Straight, full handle

Trapezoid Knives



ITEM NO.	WIDTH	DESCRIPTION
55-1214	1.2mm–1.4mm	Angled, full handle
55-1416	1.4mm–1.6mm	Angled, full handle
55-1520	1.5mm–2.0mm	Angled, full handle
55-2724	2.75mm–3.2mm	Angled, full handle
55-2824	2.85mm–3.2mm	Angled, full handle

Sideport Knives-MVR Design



ITEM NO.	WIDTH	DESCRIPTION
58-2008	0.8 mm	Angled, full handle
58-2010	1.0mm	Angled, full handle
58-2012	1.2mm	Angled, full handle
58-2014	1.4mm	Angled, full handle

Implant Knives



ITEM NO.	WIDTH	DESCRIPTION
56-5361	3.2mm	Angled, full handle, bevel up
56-3561	3.5mm	Angled, full handle, bevel up
56-4061	4.0mm	Angled, full handle, bevel up
56-5261	5.2mm	Angled, full handle, bevel up
56-5561	5.5mm	Angled, full handle, bevel up
56-6261	6.2mm	Angled, full handle, bevel up

Specialty Knives



ITEM NO.	DESCRIPTION
58-5700	Scleral, straight, blade only
58-6400	Round Tip, straight, blade only
58-6600	Lamellar, angled, blade only
58-6700	Mini, straight, blade only
58-6810	Needle blade, straight, full handle
58-6900	All Around, straight, blade only

Disc/Spoon Knives



ITEM NO.	WIDTH	DESCRIPTION
54-2261	2.2mm	Angled, full handle, bevel up
54-3061	3.0mm	Angled, full handle, bevel up

Myringotomy Blades



ITEM NO.	DESCRIPTION
58-7100	45°, downward cut, offset
58-7110	45°, upward cut, offset
58-7120	Spear Tip, narrow shaft, offset
58-7121	Lance Tip, narrow shaft, offset

To locate a UniqueEdge® Distributor or for ordering information, contact us at www.unique-tech-inc.com or by Phone at 610-775-9191

Made in the USA

- 100%
- Satisfaction guaranteed
- ISO and CE certification confirms consistency and reliability of manufacturing processes and product quality

Private label opportunities

- We manufacture for some of the largest and most trusted names in the ophthalmic industry
- Contact us with your specific product needs

ISO 9001-2000 certified & ISO 13485-2003 certified

CE Mark Authorized under MDD 93/42/EEC

unique  **technologies inc.**

Mohnton PA USA

Office: 610-775-9191 Fax: 610-775-9192

 www.unique-tech-inc.com